


## M•Biopsy™ Semi-Automatic Biopsy Instrument

The M•Biopsy™ Semi-Automatic, spring loaded biopsy instrument is compact, lightweight, disposable and adjustable.

An optional “sure-fit” coaxial introducer needle is available packaged separately (N302) or packaged in combination with the matching M•Biopsy™ Semi-Automatic Biopsy Instrument (NS301).

### Features:

- **Sharp cannula** for safe and painless percutaneous insertion and penetration.
- **Color coded** for easy gauge size determination.
- **Adjustable specimen notch** adjusts exposing either 10mm or 20mm slot providing clinical flexibility.
- **“Sure-fit” design.**
- **Echogenic cannula tip** for correct positioning under ultrasound guidance.
- **Centimeter markings** for optimal penetration control.
- **Specimen notch indicator** on handle for easy identification



Specimen notch fully exposed



Specimen notch closed



Coaxial introducer needle in combination with the biopsy instrument



Full Matrix on page 2

# M°Biopsy™ Semi-Automatic Biopsy Instrument



| M°Biopsy™                        |       |        |         |
|----------------------------------|-------|--------|---------|
| Semi-Automatic Biopsy Instrument |       |        |         |
| Item no.                         | Gauge | Length | Qty/box |
| N30114006                        | 14    | 6 CM   | 10      |
| N30114008                        | 14    | 8 CM   | 10      |
| N30114010                        | 14    | 10 CM  | 10      |
| N30114012                        | 14    | 12 CM  | 10      |
| N30114016                        | 14    | 16 CM  | 10      |
| N30114020                        | 14    | 20 CM  | 10      |
| N30114025                        | 14    | 25 CM  | 10      |
| N30116006                        | 16    | 6 CM   | 10      |
| N30116008                        | 16    | 8 CM   | 10      |
| N30116010                        | 16    | 10 CM  | 10      |
| N30116012                        | 16    | 12 CM  | 10      |
| N30116016                        | 16    | 16 CM  | 10      |
| N30116020                        | 16    | 20 CM  | 10      |
| N30116025                        | 16    | 20 CM  | 10      |
| N30118006                        | 18    | 6 CM   | 10      |
| N30118008                        | 18    | 8 CM   | 10      |
| N30118010                        | 18    | 10 CM  | 10      |
| N30118012                        | 18    | 12 CM  | 10      |
| N30118016                        | 18    | 16 CM  | 10      |
| N30118020                        | 18    | 20 CM  | 10      |
| N30118025                        | 18    | 25 CM  | 10      |
| N30118030                        | 18    | 30 CM  | 10      |
| N30120006                        | 20    | 6 CM   | 10      |
| N30120008                        | 20    | 8 CM   | 10      |
| N30120010                        | 20    | 10 CM  | 10      |
| N30120012                        | 20    | 12 CM  | 10      |
| N30120016                        | 20    | 16 CM  | 10      |
| N30120020                        | 20    | 20 CM  | 10      |

| M°Biopsy™      |       |         |         |
|----------------|-------|---------|---------|
| Coaxial Needle |       |         |         |
| Item no.       | Gauge | Length  | Qty/box |
| N30214008      | 13    | 3.2 CM  | 10      |
| N30214010      | 13    | 5.2 CM  | 10      |
| N30214016      | 13    | 11.2 CM | 10      |
| N30214020      | 13    | 15.2 CM | 10      |
| N30216008      | 15    | 3.2 CM  | 10      |
| N30216010      | 15    | 5.2 CM  | 10      |
| N30216016      | 15    | 11.2 CM | 10      |
| N30216020      | 15    | 15.2 CM | 10      |
| N30218008      | 17    | 3.2 CM  | 10      |
| N30218010      | 17    | 5.2 CM  | 10      |
| N30218012      | 17    | 7.2 CM  | 10      |
| N30218016      | 17    | 11.2 CM | 10      |
| N30218020      | 17    | 15.2 CM | 10      |
| N30220008      | 19    | 3.2 CM  | 10      |
| N30220010      | 19    | 5.2 CM  | 10      |
| N30220012      | 19    | 7.2 CM  | 10      |
| N30220016      | 19    | 11.2 CM | 10      |
| N30220020      | 19    | 15.2 CM | 10      |

| M°Biopsy™  |       |        |         |
|--|-------|--------|---------|
| Semi-Automatic Biopsy Instrument with Coaxial Needle |       |        |         |
| Item no.   | Gauge | Length | Qty/box |
| NS3011410  | 14    | 10 CM  | 10      |
| NS3011416  | 14    | 16 CM  | 10      |
| NS3011420  | 14    | 20 CM  | 10      |
| NS3011610  | 16    | 10 CM  | 10      |
| NS3011616  | 16    | 16 CM  | 10      |
| NS3011620  | 16    | 20 CM  | 10      |
| NS3011810  | 18    | 10 CM  | 10      |
| NS3011816  | 18    | 16 CM  | 10      |
| NS3011820  | 18    | 20 CM  | 10      |
| NS3012010  | 20    | 10 CM  | 10      |
| NS3012016  | 20    | 16 CM  | 10      |
| NS3012020  | 20    | 20 CM  | 10      |

