



PREMIUM AORTIC PUNCHES

Providing surgeons with the precise performance they demand during CABG surgery

TELEFLEX PREMIUM DIAMOND EDGE AORTIC PUNCH

FREE-FLOAT CUTTING ACTION TECHNOLOGY

By allowing the diamond edge blade to “float” around the punch core during operation, alignment is maintained and mechanical jamming is prevented. One punch can be used to create multiple punches in a single procedure.

DUAL-CUTTING ACTION

Two diamond-sharp cutting surfaces lift and cleanly cut through all layers of tough tissue. The smooth cutting action reduces trauma while the tapered distal portion increases visibility around the anastomotic site.

PILLING® DISPOSABLE AORTIC PUNCH – THE CLEAR ADVANTAGE®

PRECISE, CLEAN CUTS

The anvil design ensures easy insertion, while the smooth action shaft and cast stainless steel punch tips provide clean and accurate cuts, even through the most calcified tissue.

UNBEATABLE CONTROL

The countersunk thumb seat on the plunger reduces the risk of slippage, and control is enhanced by the use of specially designed finger grips. Precision tolerances minimise the risk of jamming.



PREMIUM DIAMOND-EDGE AORTIC PUNCH, LATEX FREE

Teleflex


REF.			
FULL LENGTH STERILE (7" - 3/4" HANDLE)		MID-LENGTH STERILE (6" - 3/4" HANDLE)	
		SIZE	QTY
DP-28K	MDP-28K	2.8 mm	6
DP-36K	MDP-36K	3.6 mm	6
DP-40K	MDP-40K	4.0 mm	6
DP-44K	MDP-44K	4.4 mm	6
DP-48K	MDP-48K	4.8 mm	6
DP-52K	MDP-52K	5.2 mm	6
DP-56K	MDP-56K	5.6 mm	6
DP-60K	MDP-60K	6.0 mm	6



CLEAR ADVANTAGE DISPOSABLE AORTIC PUNCH, LATEX FREE

Pilling

REF.		
FULL LENGTH STERILE (7" HANDLE) CODE NUMBER	SIZE	QTY
353427	2.7 mm	6
353435	3.5 mm	6
353440	4.0 mm	6
353445	4.5 mm	6
353447	4.7 mm	6
353450	5.0 mm	6
353455	5.5 mm	6



DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND · IDA Business Park · Athlone · Co Westmeath

Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflexmedical.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflexmedical.com

SOUTH AFRICA Tel. +27 (0)11 807 4887 · assist.africa@teleflexmedical.com

WWW.TELEFLEXMEDICAL.COM

94 05 78 · 00 00 01 · REV A · 07 10 02 · All data current at time of printing. Subject to technical changes without further notice.

Teleflex
MEDICAL