


Corneal Marking Pad

Aspen Surgical's Corneal Marking Pad can be used with many refractive/corneal marking instruments. This pre-inked, ready-to-use gentian violet marking pad saves prep time and allows you to focus on your patient. Sterile, single-use, and latex-free.



Part No.	Description	Quantity Per Box
40490	Corneal Marking Pad	10

Dr. Fog® Original Anti-Fog

DF-3100, the original anti-fog solution, comes with our easy squeeze bottle and sponge kit. We also offer the DF-3120 Ready-To-Use Sponge which is dry to touch. The sponge is alcohol-free. Both configurations are sterile and latex-free.



Part No.	Description	Quantity Per Box
DF-3100	Original Solution with Sponge	20
DF-3120	Alcohol-Free Dry Treated Sponge	40

Endoscopic Kittner


Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft. It is available in both single and dual tip. The radiopaque tip is non-absorbent. The brushed finish allows for safe use during laser surgery. Sterile, single-use, and latex-free.



Part No.	Description	Quantity Per Box
KT-9101	Single Tip	20
KT-9105	Dual Tip	20

Cautery Tip Cleaner

Cautery Tip Cleaners assist in the cleaning and removing of material from cautery tips during surgical procedures. Our tip cleaner consists of a 2" x 2" foam pad with abrasive material on one side and strong adhesive on the other for attachment to both disposable and non-disposable drapes. It will clean both mono-polar and bi-polar probes, is radiopaque, latex-free, and intended for single-use.



Part No.	Quantity Per Box
032001pbx	35
TC-100	100
TC-100b	240