


Vessel Loops

Vessel loops are used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessel loops are extruded from medical-grade silicone, with the exception of the Richard-Allan vessel loops, which are made of high performance, medical grade elastomer. All vessel loops are latex-free, radiopaque, and intended for single-use.



2 Loops / Pkg; 10 Pkgs / Box Double Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
011001pbx	Red Mini	0.83 x 406	0.03 x 16
011002pbx	Blue Mini	0.83 x 406	0.03 x 16
011003pbx	Yellow Mini	0.83 x 406	0.03 x 16
011004pbx	White Mini	0.83 x 406	0.03 x 16
011011pbx	Red Maxi	1.0 x 406	0.04 x 16
011012pbx	Blue Maxi	1.0 x 406	0.04 x 16
011013pbx	Yellow Maxi	1.0 x 406	0.04 x 16
011014pbx	White Maxi	1.0 x 406	0.04 x 16
011021pbx	White Supermaxi	1.3 x 559	0.05 x 22
011022pbx	Blue Supermaxi	1.3 x 559	0.05 x 22
011032pbx	Blue Micro	0.41 x 559	0.019 x 22
011301pbx	Red Mini	0.83 x 762	0.03 x 30
011302pbx	Blue Mini	0.83 x 762	0.03 x 30
011303pbx	Yellow Mini	0.83 x 762	0.03 x 30
011304pbx	White Mini	0.83 x 762	0.03 x 30
011311pbx	Red Maxi	1.0 x 762	0.04 x 30
011312pbx	Blue Maxi	1.0 x 762	0.04 x 30
011313pbx	Yellow Maxi	1.0 x 762	0.04 x 30
011314pbx	White Maxi	1.0 x 762	0.04 x 30

Vessel Loops continued

2 Loops / Pkg; 10 Pkgs / Box Carded and Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
016001pbx	Red Mini	0.83 x 457	0.03 x 18
016002pbx	Blue Mini	0.83 x 457	0.03 x 18
016003pbx	Yellow Mini	0.83 x 457	0.03 x 18
016004pbx	White Mini	0.83 x 457	0.03 x 18
016011pbx	Red Maxi	1.0 x 457	0.04 x 18
016012pbx	Blue Maxi	1.0 x 457	0.04 x 18
016013pbx	Yellow Maxi	1.0 x 457	0.04 x 18
016014pbx	White Maxi	1.0 x 457	0.04 x 18

RICHARD-ALLAN Vessel Loops



2 Loops / Pkg; 50 Pkgs / Box Banded and Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
3901	Maxi Red	1.0 x 457	0.04 x 18
3902	Maxi Blue	1.0 x 457	0.04 x 18
3903	Maxi White	1.0 x 457	0.04 x 18
3904	Maxi Yellow	1.0 x 457	0.04 x 18
3905	Mini Red	.83 x 457	0.03 x 18
3906	Mini Blue	.83 x 457	0.03 x 18
3907	Mini White	.83 x 457	0.03 x 18
3908	Mini Yellow	.83 x 457	0.03 x 18

4 Loops / Pkg; 25 Pkgs / Box

Part No.	Color	mm Width x Length	Inches Width x Length
3910	Four Color Maxi Set	1.0 x 457	0.04 x 18
3912	Four Color Mini Set	0.83 x 457	0.03 x 18